

- PanCO: Results of a single-arm pilot study of <sup>32</sup>P microparticles in unresectable locally advanced pancreatic adenocarcinoma
   with gemcitabine/nab-paclitaxel or FOLFIRINOX chemotherapy<sup>1</sup>
  - P J Ross, H S Wasan, D Croagh, M Nikfarjam, N Nguyen, M Aghmesheh, A M Nagrial, D Bartholomeusz, A Hendlisz, T Ajithkumar, C Iwuji, N E Wilson, D M Turner, D C James, E Young, M T Harris



## **Study Objective**

To further investigate the safety, efficacy, feasibility and performance of the OncoSil™ device when implanted intratumourally using EUS in a patient population undergoing standard chemotherapy for unresectable LAPC.

## **Study Design**



#### **Key Eligibility Criteria**

- Histologically or cytologically proven pancreatic adenocarcinoma
- Unresectable LAPC
- Target tumour diameter 2-6cm
- ECOG Performance Status 0 to 1
- · No distant metastases
- No prior radiotherapy or chemotherapy for pancreatic cancer



#### Location

International, multicentre pilot study with 10 sites in 3 countries:

- Australia
- Belgium
- UK



#### Data Collection/Follow-up

- SPECT-CT Bremsstrahlung imaging at ≤4 hours and day 7
- Blood and urine <sup>32</sup>P analysis
- 8 weekly CT RECIST 1.1 and tumour volume\*
- FDG-PET Baseline and week 12\*
- CA 19-9 tumour marker serial analysis



#### Treatment

<sup>32</sup>P activity calculated from patients' tumour volume to deliver 100 Gy absorbed dose.



WEEK 4 OncoSil™ Device Implantation WEEK 5 ONWARDS
Continuation
of
chemotherapy


PD

Overall survival

† By physician choice; per standard-of-care



## **Participant Flow**





## **Baseline Characteristics**

Demographic/Characteristic		n (%) ITT Cohort (N=50)
Age, years	Median (Range)	65 (42-84)
Sex	Male : Female	28 (56%) : 22 (44%)
Race:	White/Caucasian Asian Black/African American	40 (80%) 7 (14%) 3 (6%)
ECOG Performance Status 0:1		26 (52%) : 24 (48%)
CA 19-9, (U/	mL) [n=49] Median (Range)	163 (1–6576)
Pancreatic tumour location Head : Body		42 (84%) : 8 (16%)
Target lesion longest diameter, cm* Median (Range)		4.5 (2.6-7.1)
Tumour volume, cc* Median (Range)		24.4 (7.9-68.7)
Study days to OncoSil™ implantation [n=42] Median (Range) 31 (21-77)		
Chemother	gemcitabine + nab-paclitaxel FOLFIRINOX	40 (80%) 10 (20%)

\*By independent central reader analysis Oncosil.com





## Safety and Tolerability (PP Cohort)

- 956 TEAEs reported; 139 were Grade ≥3
   No serious device or radiation-related toxicities were reported
- 289 TEAEs (30.2%) occurred pre-OncoSil<sup>™</sup> implantation (25 Grade ≥3) vs. 667 TEAEs (69.8%) post-implant (114 Grade ≥3) [median follow-up: 1 vs. 31 months, respectively], with 41 vs. 609 attributed to the OncoSil<sup>™</sup> device and/or implantation procedure vs. chemotherapy, respectively
- No increased incidence of grade ≥3 TEAEs pre-implantation vs. per cycle, overall and by key categories, post-implantation

OncoSil™ in combination
with first-line chemotherapy
for LAPC had an
acceptable safety profile
over a study timeframe
(median 31.6 months)



Tumour Response (Best Response by Central Imaging Analysis, RECIST v1.1\*^)

Best Response, Evaluable Patients	ITT Cohort (n=50)	PP Cohort (n42)
Complete Response [CR]**	0 (0%)	0 (0%)
Partial Response [PR]**	14 (29.8%)	13 (31.0%)
Stable Disease [SD]**	31 (66%)	29 (69.0%)
Progressive Disease [PD]**	2 (4.3%)	0 (0%)
Not evaluated	3	0
Overall Response Rate [ORR]	14 (28.0%)	13 (31.0%)
Disease Control Rate [DCR]	45 (90.0%)	42 (100%)

<sup>\*</sup>By central image reader analysis. \*\*Percentages based on the number of assessable study participants. \*Response before surgical resection.



### Local Disease Control Rate (LDCR<sub>16 weeks</sub>) and Surgical Resection

Outcome	ITT Cohort (n=50) 41 (82.0%)	PP Cohort (n=42) 38 (90.5%)
LDCR <sub>16 weeks</sub> , n (%) [1° Efficacy Endpoint]		
Rate of Surgical Resection with Curative Intent, n (%)	10 (20.0%)	10 (23.8%)
<b>RO Margins, n (%)</b> vs. R1 Margins, n (%)	<b>8 (80.0%)</b> vs. 2 (20.0%)	<b>8 (80.0%)</b> vs. 2 (20.0%)
Median Progression-Free Survival (95% CI)	9.3 months (5.7-11.3)	9.3 months (5.8-11.3)
<b>Median Overall Survival</b> (95% CI)	<b>15.2 months</b> (11.3-18.8)	<b>15.5 months</b> (11.4-20.1)

#### Abbreviations:

CI: Confidence interval

**ECOG:** Eastern Cooperative Oncology Group

**EUS:** Endoscopic ultrasound

ITT: Intention-to-treat

LAPC: Locally advanced pancreatic cancer

NC: Not calculable

## 33% patients

(14/42) were either surgically resected OR were technically sufficiently downstaged to be considered for resection\*



# 23.8% patients (10/42) were resected

- 9 patients received gemcitabine + nab-paclitaxel; 1 received FOLFIRINOX
- HPB surgeons noted reduction in the fibrosis of the tumours along blood vessels and favourable tissue planes
- 4 resected patients died (at 18.8-19.2 months from enrolment)
  - 6 remained alive at 26.4-35.3 months post-enrolment (5 disease-free)

\*Four patients did not undergo surgery due to concomitant co-morbidities and/or other considerations (advanced age, patient choice).

**PD:** Progressive disease **PP:** Per protocol

TEAEs: Treatment-emergent adverse events

REFERENCES: 1. Ross PJ, Wasan HS, Croagh D et al. Results of a single-arm pilot study of 32° microparticles in unresectable locally advanced pancreatic adenocarcinoma with gemcitabine/nab-paclitaxel or FOLFIRINOX chemotherapy. ESMO Open February 2022;7(1):100356. doi: 10.1016/j.esmoop.2021.

2. Ross PJ, Burnett C, Nikfarjam M et al. Comparison of Resected vs. Non-resected PanCO Study Patients with Unresectable Locally Advanced Pancreatic Cancer (uLAPC) Receiving 32 P-microparticles and Chemotherapy. Poster presentation PBP-006 E-AHPBA 2023, Lyon, France. https://doi.org/10.1016/j.hpb.2023.07.454.

INTENDED USE / INDICATIONS FOR USE: OncoSil™ is intended for intratumoural implantation into a pancreatic tumour via injection under endoscopic ultrasound guidance. OncoSil™ is indicated for the treatment of patients with locally advanced unresectable pancreatic cancer, in combination with gemcitabine-based chemotherapy.

This information is intended for healthcare professionals only. All medical treatments carry benefits and risks. For safety related information, please refer to the OncoSil® System Instructions for Use.